



ERAL1 Antibody

Product Code	CSB-PA007759GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75616
Immunogen	Human ERAL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Era G-protein-like 1 (E. coli);ERAL1;ERA;ERAL1A;HERA-A;HERA-B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERAL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.