



ERAP2 Antibody

Product Code	CSB-PA007761GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P179
Immunogen	Human ERAP2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	endoplasmic reticulum aminopeptidase 2;ERAP2;FLJ23633;FLJ23701;FLJ23807;L-RAP;LRAP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.