



EXOSC10 Antibody

Product Code	CSB-PA007889GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01780
Immunogen	Human EXOSC10
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	exosome component 10;EXOSC10;PM-Sc1;PM/Sc1-100;PMSCL;PMSCL2;RRP6;Rrp6p;p2;p3;p4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EXOSC10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.