



EXOSC4 Antibody

Product Code	CSB-PA007892GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPD3
Immunogen	Human EXOSC4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	exosome component 4;EXOSC4;FLJ20591;RRP41;RRP41A;Rrp41p;SKI6;Ski6p;hRrp41p;p12A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EXOSC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.