



FA2H Antibody

Product Code	CSB-PA007937GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L5A8
Immunogen	Human FA2H
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	fatty acid 2-hydroxylase;FA2H;FAAH;FAH1;FAXDC1;FLJ25287;SCS7 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FA2H
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.