



FAM107A Antibody

Product Code	CSB-PA007987GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95990
Immunogen	Human FAM107A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	family with sequence similarity 107, member A;FAM107A;DRR1;FLJ30158;FLJ45473;TU3A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FAM107A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.