



FAM117B Antibody

Product Code	CSB-PA008019GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1L5
Immunogen	Human FAM117B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	family with sequence similarity 117, member B;FAM117B;ALS2CR13;DKFZp686H01244;FLJ38771 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FAM117B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.