



FAM175A Antibody

Product Code	CSB-PA008138GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UWZ7
Immunogen	Human FAM175A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	family with sequence similarity 175, member A;FAM175A;ABRA1;CCDC98;FLJ11520;FLJ12642;FLJ13614 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABRAXAS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.