



HAUS6 Antibody

Product Code	CSB-PA008208GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z4H7
Immunogen	Human FAM29A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	family with sequence similarity 29, member A; FAM29A; Dgt6; MGC102696; MGC138798; MGC138799; RP11-296P7.3 ; DGT6, HAUS6, KIAA1574
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HAUS6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.