



RIPOR2 Antibody

Product Code	CSB-PA008285GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4F9
Immunogen	Human FAM65B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	family with sequence similarity 65, member B;FAM65B;C6orf32;DIFF40;DIFF48;KIAA0386;PL48 ;
Immunogen Species	Homo sapiens (Human)
Target Names	RIPOR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.