



FBF1 Antibody

Product Code	CSB-PA008449GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TES7
Immunogen	Human FBF1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Fas (TNFRSF6) binding factor 1;FBF1;Alb;FBF-1;FLJ00103;FLJ13996;KIAA1863 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FBF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.