



FCER1A Antibody

Product Code	CSB-PA008532GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12319
Immunogen	Human FCER1A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Fc fragment of IgE, high affinity I, receptor for; alpha polypeptide;FCER1A;FCE1A;FcERI ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCER1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.