



FCRLB Antibody

Product Code	CSB-PA008561GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6BAA4
Immunogen	Human FCRLB
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Fc receptor-like B;FCRLB;FCRL2;FCRLM2;FCRLY;FLJ31052;FREB-2;FcRY;RP11-474I16.6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCRLB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.