



# FERMT2 Antibody

<b>Product Code</b>	CSB-PA008594GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96AC1
<b>Immunogen</b>	Human FERMT2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	fermitin family homolog 2 (Drosophila);FERMT2;DKFZp686G11125;FLJ34213;FLJ44462;KIND2;MIG2;PLEKHC1;UNC112;UNC112B;mig-2 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FERMT2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.