



FKBP14 Antibody

Product Code	CSB-PA008689GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWM8
Immunogen	Human FKBP14
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	FK506 binding protein 14, 22 kDa;FKBP14;FKBP22;FLJ20731 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FKBP14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.