



# FMO3 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA008749GA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                        |
| <b>Uniprot No.</b>         | P31513   |
| <b>Immunogen</b>           | Human FMO3   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA,WB,IHC   |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Purification Method</b> | Ammonium Sulfate Precipitation   |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | flavin containing monooxygenase 3;FMO3;FMOII;MGC34400;TMAU;dJ127D3.1<br>;            |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                       |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | FMO3   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.          |