



FOXM1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA008828GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q08050 |
| Immunogen | Human FOXM1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | forkhead box M1;FOXM1;FKHL16;FOXM1B;HFH-11;HFH11;HNF-3;INS-1;MPHOSPH2;MPP-2;MPP2;PIG29;TGT3;TRIDENT ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FOXM1 |