



FOXM1 Antibody

Product Code	CSB-PA008828GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08050
Immunogen	Human FOXM1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	forkhead box M1;FOXM1;FKHL16;FOXM1B;HFH-11;HFH11;HNF-3;INS-1;MPHOSPH2;MPP-2;MPP2;PIG29;TGT3;TRIDENT ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.