



FRS3 Antibody

Product Code	CSB-PA009005GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43559
Immunogen	Human FRS3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	fibroblast growth factor receptor substrate 3;FRS3;FRS2B;FRS2beta;MGC17167;SNT-2;SNT2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FRS3
Target Details	This protein is a substrate for the fibroblast growth factor receptor. It is found in peripheral plasma membrane and functions in linking FGF receptor stimulation to activators of Ras.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.