



FSD1L Antibody

Product Code	CSB-PA009014GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXM9
Immunogen	Human FSD1L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	fibronectin type III and SPRY domain containing 1-like;FSD1L;CCDC10;CSDUFD1;FSD1CL;FSD1NL;MGC45564;MIR1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FSD1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.