



GABBR2 Antibody

Product Code	CSB-PA009135GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75899
Immunogen	Human GABBR2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	gamma-aminobutyric acid (GABA) B receptor, 2;GABBR2;FLJ36928;GABABR2;GPR51;GPRC3B;HG20;HRIHFB2099 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GABBR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.