



GADD45GIP1 Antibody

Product Code	CSB-PA009165GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAE8
Immunogen	Human GADD45GIP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	growth arrest and DNA-damage-inducible, gamma interacting protein 1;GADD45GIP1;CKBBP2;CRIF1;MGC4667;MGC4758;PLINP-1;PRG6;Plinp1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GADD45GIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.