



GBP5 Antibody

Product Code	CSB-PA009301GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PP8
Immunogen	Human GBP5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	guanylate binding protein 5;GBP5;GBP-5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GBP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.