



GHRL Antibody

Product Code	CSB-PA009414GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBU3
Immunogen	Human GHRL
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ghrelin/obestatin prepropeptide;GHRL;MTLRP;obestatin ;
Immunogen Species	Homo sapiens (Human)
Target Names	GHRL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.