



GIPC1 Antibody

Product Code	CSB-PA009435GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14908
Immunogen	Human GIPC1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	GIPC PDZ domain containing family, member 1;GIPC1;C19orf3;GIPC;GLUT1CBP;Hs.6454;IIP-1;MGC15889;MGC3774;NIP;RGS19IP1;SEMCAP;SYNECTIIN;SYNECTIN;TIP-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GIPC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.