



GLUL Antibody

Product Code	CSB-PA009553GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15104
Immunogen	Human GLUL
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	glutamate-ammonia ligase (glutamine synthetase);GLUL;GLNS;GS;PIG43;PIG59 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GLUL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.