



GLYATL1 Antibody

Product Code	CSB-PA009559GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969I3
Immunogen	Human GLYATL1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	glycine-N-acyltransferase-like 1; GLYATL1; FLJ26507; FLJ34646; GATF-C; GNAT; MGC15397; MGC15937 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GLYATL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.