



GNAQ Antibody

Product Code	CSB-PA009594GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50148
Immunogen	Human GNAQ
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	guanine nucleotide binding protein (G protein), q polypeptide;GNAQ;G-ALPHA-q;GAQ ;
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNAQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.