



GNB2L1 Antibody

Product Code	CSB-PA009605GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63244
Immunogen	Human GNB2L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	guanine nucleotide binding protein (G protein), beta polypeptide 2-like 1;GNB2L1;Gnb2-rs1;H12.3;HLC-7;PIG21;RACK1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RACK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.