



GNL3 Antibody

Product Code	CSB-PA009625GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BVP2
Immunogen	Human GNL3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	guanine nucleotide binding protein-like 3 (nucleolar); GNL3; C77032; E2IG3; MGC800; NS ; Nucleostemin; NNP47
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.