



GPR132 Antibody

Product Code	CSB-PA009749GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNW8
Immunogen	Human GPR132
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	G Protein Coupled Receptor G2A antibody; G2 accumulation protein antibody; G2A antibody; GP132_HUMAN antibody; GPR132 antibody; Probable G-protein coupled receptor 132 antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPR132
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.