



Recombinant Human Eukaryotic initiation factor 4A-III (EIF4A3)

| | |
|--------------------------|--|
| Product Code | CSB-MP007553HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P38919 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATTATMATS GSARKRLLKE EDMTKVEFET SEEVDVTPTF DTMGLREDLL RGIYAYGF EK PS AIQQRAIK QIIKGRDVIA QSQSGTGKTA TFSISVLQCL DIQVRETQAL ILAPTRELAV QIQKGLLALG DYMNVQCHAC IGGTNVGEDI RKLDYGQHV V AGTPGRV FDM IRRRSLRTRA IKMLVLDEAD EMLNKGFKEQ YDVYRYLPP ATQVVLISAT LPHEILEMTN KFMTDPIRIL VKRDEL TLEG IKQFFVAVER EEWKFD TLCD LYDTLTITQA VIFCNTKRKV DWLTEKMREA NFTVSSMHGD MPQKERESIM KEFRSGASRV LISTDVWARG LDVPQVSLII NYDLPNNREL YIHRIGRSGR YGRKGVAINF VKNDDIRILR DIEQYYSTQI DEMPMNVADL I |
| Source | Mammalian cell |
| Target Names | EIF4A3 |
| Protein Names | Recommended name: Eukaryotic initiation factor 4A-III Short name= eIF-4A-III Short name= eIF4A-III EC= 3.6.4.13 Alternative name(s): ATP-dependent RNA helicase DDX48 ATP-dependent RNA helicase eIF4A-3 DEAD box protein |
| Expression Region | 1-411 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This gene encodes a member of the DEAD box protein family. DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure, such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This protein is a nuclear matrix protein. Its amino acid sequence is highly similar to the amino acid sequences of the translation initiation factors eIF4AII and eIF4AIII, two other members of the DEAD box protein family. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.