



GTF2E2 Antibody

Product Code	CSB-PA010003GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29084
Immunogen	Human GTF2E2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	general transcription factor IIE, polypeptide 2, beta 34kDa;GTF2E2;FE;TF2E2;TFIIE-B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Transcription
Target Names	GTF2E2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.