



GTF3C4 Antibody

Product Code	CSB-PA010028GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKN8
Immunogen	Human GTF3C4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	general transcription factor IIIC, polypeptide 4, 90kDa;GTF3C4;FLJ21002;KAT12;MGC138450;TFIII90;TFIIIC90;TFIIICdelta;TFi iiC2-90 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTF3C4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.