



GTPBP3 Antibody

Product Code	CSB-PA010034GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969Y2
Immunogen	Human GTPBP3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	GTP binding protein 3 (mitochondrial);GTPBP3;FLJ14700;GTPBG3;MSS1;MTGP1;THDF1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GTPBP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.