



MTG2 Antibody

Product Code	CSB-PA010036GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4K7
Immunogen	Human GTPBP5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	GTP binding protein 5 (putative);GTPBP5;FLJ10741;MGC29512;ObgH1;dJ1005F21.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MTG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.