



HPD Antibody

Product Code	CSB-PA010698GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32754
Immunogen	Human HPD
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	4-hydroxyphenylpyruvate dioxygenase;HPD;4-HPPD;4HPPD;GLOD3;PPD ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HPD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.