



ICA1L Antibody

Product Code	CSB-PA010948GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDH6
Immunogen	Human ICA1L
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	islet cell autoantigen 1,69kDa-like; ICA1L; ALS2CR14; ALS2CR15; DKFZp434E1919; MGC138440 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ICA1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.