



ICOSLG Antibody

Product Code	CSB-PA010958GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75144
Immunogen	Human ICOSLG
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	inducible T-cell co-stimulator ligand; ICOSLG; B7-H2; B7H2; B7RP-1; B7RP1; CD275; GL50; ICOS-L; ICOSL; KIAA0653; LICOS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ICOSLG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.