



IFIT5 Antibody

Product Code	CSB-PA011023GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13325
Immunogen	Human IFIT5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	interferon-induced protein with tetratricopeptide repeats 5; IFIT5; FLJ53857; FLJ92678; RI58 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFIT5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.