



IGF2BP3 Antibody

Product Code	CSB-PA011092GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00425
Immunogen	Human IGF2BP3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	insulin-like growth factor 2 mRNA binding protein 3;IGF2BP3;DKFZp686F1078;IMP-3;IMP3;KOC1;VICKZ3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IGF2BP3
Target Details	This protein is primarily found in the nucleolus, where it can bind to the 5 UTR of the insulin-like growth factor II leader 3 mRNA and may repress translation of insulin-like growth factor II during late development. The encoded protein contains several KH domains, which are important in RNA binding and are known to be involved in RNA synthesis and metabolism. A pseudogene exists on chromosome 7, and there are putative pseudogenes on other chromosomes.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.