



IGKV1-5 Antibody

Product Code	CSB-PA011359GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01602
Immunogen	Human IGKV1-5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	immunoglobulin kappa variable 1-5;IGKV1-5;IGKV;IGKV15;L12;L12a;MGC22745;MGC32715;MGC88810;V1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IGKV1-5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.