



IKBKG Antibody

Product Code	CSB-PA011574GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K9
Immunogen	Human IKBKG
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase gamma;IKBKG;AMCBX1;FIP-3;FIP3;Fip3p;IKK-gamma;IP;IP1;IP2;IPD2;NEMO ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IKBKG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.