



# IL11RA Antibody

<b>Product Code</b>	CSB-PA011584GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q14626
<b>Immunogen</b>	Human IL11RA
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Storage Buffer</b>	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IL11RA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.