



# IL17A Antibody

<b>Product Code</b>	CSB-PA011597GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q16552
<b>Immunogen</b>	Human IL17A
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,IHC
<b>Form</b>	Liquid
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	interleukin 17A;IL17A;CTLA8;IL-17;IL-17A;IL17 ;
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IL17A
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.