



IL17RB Antibody

Product Code	CSB-PA011603GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRM6
Immunogen	Human IL17RB
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	interleukin 17 receptor B; IL17RB; CRL4; EVI27; IL17BR; IL17RH1; MGC5245 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL17RB
Target Details	This protein is a cytokine receptor. This receptor specifically binds to IL17B and IL17E, but does not bind to IL17 and IL17C. This receptor has been shown to mediate the activation of NF-kappaB and the production of IL8 induced by IL17E. The expression of the rat counterpart of this gene was found to be significantly up-regulated during intestinal inflammation, which suggested the immunoregulatory activity of this receptor.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.