



IL17RE Antibody

Product Code	CSB-PA011606GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFR9
Immunogen	Human IL17RE
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	interleukin 17 receptor E;IL17RE;FLJ23658;MGC71884 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL17RE
Target Details	This gene encodes a protein that shares limited similarity with the interleukin-17 receptor. The function of the encoded protein has not yet been determined. Three alternatively spliced transcript variants of this gene encoding two distinct isoforms have been reported.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.