



IL2RG Antibody

Product Code	CSB-PA011651GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31785
Immunogen	Human IL2RG
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL2RG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.