



# IMP3 Antibody

<b>Product Code</b>	CSB-PA011692GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9NV31
<b>Immunogen</b>	Human IMP3
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	IMP3, U3 small nucleolar ribonucleoprotein, homolog (yeast);IMP3;BRMS2;C15orf12;DKFZp586L0118;FLJ10968;MRPS4 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IMP3
<b>Target Details</b>	This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.