



INE1 Antibody

Product Code	CSB-PA011709GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15225
Immunogen	Human INE1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	inactivation escape 1 (non-protein coding); INE1; MGC97326; MGC97335; MGC97386; NCRNA00010 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	INE1
Target Details	X chromosome inactivation provides dosage compensation for the expression level of X-linked genes from the single X in males and the two in females. This X chromosome gene was identified because its transcription escapes X inactivation in females.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.